



Realizing Successful Business Activity



NEXOBIO Co., Ltd.

The Next Generation
of Dental Solutions.

Peace for Teeth.

www.nexobio.com

PRODUCT LIST

T com	03	MTA Cem S	07
T Com Flow	03	T Seal	07
T Core	03	MTA CEM LC	08
Transeal F	04	T Pex	08
T Etchant	04	T Paste	08
T Bond	04	T EDTA Cream	09
T Block	05	Soft- Dam	09
T Hole-in	05	T Cem Implant	09
Trantemp	05	T Cem NE	10
Transcen Temp C&B LC	06	T Ortho Bond	10
Transcen Temp In&Onlay LC	06	T Ortho Band	10
T stat	06	T Com Bulk Flow	11
MTA Cem	07	T Bioseal	11

CONTENTS

PRODUCT LIST



T com

Nano-Hybrid composite Resin

- Low shrinkage stress/ Non Sticky, Non-Slumping.
- Package : • T com Kit
 - 4g syringes with 5 shade(A1, A2, A3, A3.5, B2)
 - 1 Bottle/5g of Bonding agent (T Bond)
 - 1 syringe/5g of Etching agent (T Etchant)
 - T com Refill 4g single syringe.
- Shade select : A1, A2, A3, A3.5, B1, B2, 0A2, WT



T Com Flow

Flowable Nano-Hybrid composite Resin

- 40 Nano Size Filler for excellent polishability.
- Package : • T com Flow Kit
 - 2g Syringes with 6 shade(A1, A2, A3, A3.5, B2, A2E)
 - T com Flow refill
 - 2g x 2 syringes with 10 tips.
- Shade Select : A1, A2, A3, A3.5, B1, B2, WT, A2E



T Core

Core Build-up And post Cementation

- Excellent Compressive and Flexural Strength.
- Dual Syringe type.
- Package : 9g x 2 syringes with 10 mixing tips.
- Shade Select : A3, White



Transeal F

Light Curing type Pit and Fissure Sealant

- Fluoride contained and releases.
- Good bonding strength and Low viscosity.
- Package : 1.2ml x 4 syringes.



T Etchant

37% Phosphoric Acid Etching Gel

- Conditioning of Surface for dentine and enamel
- Increase bonding strength
- Perfect rinse out completely
- Self limiting depth to only 1.9micrometers
- Package : 5g x 3 syringes with tips.



T Bond

Universal Adhesive-One Step bonding system

- Total-Etching, Selective-Etching and Self-Etching.
- Combine Wet bonding and Self adhesive.
- MDP and GPDM based coupling monomer
- Package : 5g with Micro-brushes.



T Block

Flowable Light Cured Temporary and block out

- Easy to use and remove.
- Package : 5g x 3 syringes with Applicable tips.



T Hole-in

Light Curing Semi-gel abutment Hole filling

- Implant-Screw hole filling
- Temporary filling for Inlay & Onlay
- Package : 5g x 2 syringes



Trantemp

Esthetic Temporization Materials

- Automix bis-acrylic temporary material for crowns, bridges
- Strong Temporary with Long-term Durability.
- Highest Fracture resistance.
- Package : 65g Cartridge
- Shade Select : A1, A2



Transcen Temp C&B LC

Light cured Temporary Crown& Bridge

- Fast and easy
- MMA Free
- Low Shrinkage.
- Strong & Durable
- Package: 15g x 3 sticks
- Shade Select : A2, A3



Transcen Temp In&Onlay LC

Light curing type Temporary filling

- Good adhesion and Easy remove.
- Fast Curing and High Curing depth
- Package : 2.5g x 3 syringes.
- Shade Select : Blue, White



T stat

Paste Retraction and Hemostasis System

- 15% Aluminum Chloride Topical Gingival Retraction Paste.
- Quick effective retraction- 2 minutes
- Excellence safety
- Package : 0.9g x 2 syringes with applicable tips.



MTA Cem

Mineral Trioxide Aggregate Endodontic Repair Material

- Root Repair and Pulp Capping
- Apexication and Apicotomy.
- Tooth Coloured formula and High Radiopacity.
- Package : 1g bottle / 0.14g x 3 capsules.



MTA Cem S

Root Canal Filling materials

- Excellent Biocompatibility.
- Calcium silicate ceramic based
- Package : 2.2g x 1 Syringe



T Seal

Epoxy Resin Base Root Canal Sealer

- Excellent Biocompatibility.
- Good sealing ability.
- Calcium silicate ceramic based
- 1:1 Dual Syringe type.
- Package : 12g, Spatula and Mixing Pad.



MTA CEM LC

Light Cured Calcium Silicate Pulp Protectant
Cavity Liner

- Pulp capping and Liner
- Calcium Release
- Alkaline pH-Promotes healing
- Package : 2g x 4 Syringes



T Pex

Calcium Hydroxide with Iodoform

- Weeping Canal, Apexification. Root Canal filling for deciduous teeth.
- Excellent Antibacterial effect and Radiopacity.
- Premixed Paste
- Package : 2.2g x 1 Syringe with Applicable tips.



T Paste

Calcium Hydroxide with Barium Sulfate

- Water soluble
- Easy Cleaning and removing.
- Package : 2.2g x 1 syringe with Applicable tips.



T EDTA Cream

19% EDTA, Chelating and Lubrication

- Root Canal Preparation with files and reamers.
- Smear Layer removal
- Package : 7g x 2 syringes with application tips.



Soft- Dam (Gingival Barrier)

- Protect Soft Tissue with impervious seal.
- Easy to remove.
- Minimum heat while curing.
- Package : 1.2g x 4 Syringes with Applicable tips.



T Cem Implant

Non Eugenol Temporary Cement for Implant
Luting

- Long term Cementation for implant-retained restoration and provisional.
- Setting time : 5 mins.
- Compressive resistance: 34 Mpa
- Film Thickness : 17 micro
- Package : 5ml x 2 dual syringes with mixing tips.



T Cem NE

None Euginol Temporary Cement

- mix and minimum loss.
- Setting time: 5 mins
- Package : 10g x 2 dual syringes with Mixing pad and Spatula



T Ortho Bond

Light cure Orthodontic Adhesive

- Strong reliable bonds.
- Package : 5g x 2 Syringes with Primer
5g x 1 Syringe.



T Ortho Band

- Composite resin for orthodontic band
- Package : 5g x 1 Syringe.



T Com Bulk Flow

Flowable Bulk Fill Restorative

- Maximun 4 mm increments
- self-leveling phenomenon
- Package : 2g x 3 syringes with 10 tips.
- Shade Select : A1, A2, A3.



T Bioseal

Bioceramic Root Canal Sealer

- Calcium silicate Bioceramic root canal sealing materials
- Safe and biocompatible
- Bioactivity of bioceramic sealers
- Dimensionally stable
- Package : 2g x 1 syringe with 10 tips.



NEXOBIO Co., Ltd.

1107, T1 Tower, 31, Bongmyeong-ro, Heungdeok-gu, Cheongju-si, Chungcheongbuk-do, Korea.

Tel. +82-43-298-9775 Fax. +82-43-298-9777 E-mail. nexobio@naver.com / info@nexobio.com